

## Caregiver competency checklist

Patient name:	_ Caregiver name:		
Goals		Date goals met	Trainer
The caregiver will be able to:			
Understand basic indication for AIRVO usage.			
2. Review safety Information.			
<ul> <li>3. Explain AIRVO system components.</li> <li>AIRVO</li> <li>Heated breathing tube</li> <li>Water chamber</li> <li>Patient interface</li> </ul>			
<ul> <li>4. Explain the buttons, indicators, and connectors on the frear, and side of the AIRVO.</li> <li>Display</li> <li>Power Button</li> <li>Power Cord</li> <li>Alarm Silence Button</li> <li>Heater Plate</li> <li>UP/DOWN Button</li> <li>Heated Breathing Tume Mode Button</li> <li>Oxygen supply connection</li> </ul>			
5. Attach O2 source and discuss O2 safety.			
6. Assemble and connect:  ☐ Heated breathing tube ☐ Oxygen supply hose ☐ Patient interface ☐ Water chamber			
7. Turn AIRVO on and off.			

New Jersey Respiratory Associates, Inc. 333 Highway 46, West | Suite 201-B

Ph: (973) 244-2190 Fax: (973) 244-2195 **SeeNJRA.com** 



Goals		Date goals met	Trainer
The caregiver will be able to:			
<ul> <li>8. Navigate the menu screens.</li> <li>Configure target settings:</li> <li>Dew-Point Temperature</li> <li>Target Flow</li> <li>Oxygen</li> </ul>			
9. Locate and explain alarm messages and alarm indi (audio and visual). Explain the audio pause feature Explain what to do when an alarm occurs			
10. Review care, cleaning, and maintenance of AIRVO accessories.	and		
11. Review emergency procedures.			
Caregiver signature	Technician signature		

New Jersey Respiratory Associates, Inc. 333 Highway 46, West | Suite 201-B Ph: (973) 244-2190

Fax: (973) 244-2195 SeeNJRA.com

