



# AIRVO *by Fisher & Paykel*

The AIRVO is a humidifier with integrated flow generator that delivers warmed and humidified respiratory gases to spontaneously breathing patients. NJRA provides post-acute providers everything they need to care for an individual requiring an AIRVO.

## Indications for Use

The AIRVO is for the treatment of spontaneously breathing patients who would benefit from receiving high-flow warmed and humidified respiratory gases, including patients who have upper airways bypassed. The oxygen flow may be from 2 - 60L/min depending on the patient interface and a facility's oxygen delivery source. A facility with piped-in oxygen can deliver up to 60L/min, while a facility using oxygen concentrators can deliver up to 10L/min.

## Patient Interface

The AIRVO can be used with a variety of patient interfaces. The Optiflow Nasal Interface covers the nostrils, leaving the mouth unobstructed. The Tracheostomy Interface connects directly to the tracheostomy tube and is secured using a lanyard around the neck. Finally, the Mask Interface Adapter allows users to connect the tubing to a 22mm port on a vented tracheostomy mask or standard oxygen mask. NJRA's respiratory therapist will help select the interface the best fits the individual's needs.



**Fisher & Paykel AIRVO**  
Humidifier with integrated flow generator



Optiflow Nasal Interface



Tracheostomy Interface



Adapter w/ Oxygen Mask



Adapter w/ Tracheostomy Mask

## How We Help

NJRA works with each facility to help ensure they have all the equipment and training they need to care for an individual requiring a AIRVO. Our services include:

- AIRVO Rental
- Patient interfaces, circuits, & chambers
- Equipment set-up, clinical assessment, follow-up visits
- AIRVO in-services, policy & procedures, nurse competencies

If you are interested in treating individuals using the AIRVO contact us today!